

Konelab™ dedicated reagents from MeDia diagnostici

Convenience

Dedicated reagents from MeDia offer you choice without compromise! MeDia supply a range of high specification reagents which have been 'finely tuned' to work with your Konelab™ instrument.

Versatility

Reagents are available in all sizes of dedicated packaging which are designed to fit directly onto the instrument of your choice.

Simplicity

Barcoded reagents in MeDia dedicated packaging are available to work directly with your Konelab™ software.

BudgetControl

User defined Konelab™ reagents are available from MeDia when cost considerations are paramount. Reagents can be transferred directly into Konelab™.



Konelab™ 20, T30, T60 packaging with barcode for Konelab™ defined applications.

Cat.No.	Product description	Method	Kit size	Kit.config.	No. of tests
100370460K	ALBUMIN MR	BCG	4x60	R1	1.000
101001450K	ALKALINE PHOSPHATASE	DGKC - DEA	4x40+1x40	R1+R2	1.170
101002450K	ALT (GPT)	IFCC no PDP	4x40+1x40	R1+R2	1.600
101236160K	AMYLASE MR 1X60	CNPG3	4x40	R1	500
1013364010K	ASO	LEIA	40+10	R1+R2+STD	160
101004450K	AST (GOT)	IFCC no PDP	4x40+1x40	R1+R2	1.600
100326450K	BILIRUBIN DIRECT	DCA	4x40+2x20	R1+R2	1.110
100301450K	BILIRUBIN TOTAL	DCA	4x40+2x20	R1+R2	1.110
100272460K	CALCIUM ASX	Arsenazo III	4x60	R1	1.400
101295460K	CHOLESTEROL TOTAL	CHOD-PAP	4x60	R1	1.090
101006520K	CK NAC activated	DGKC	5x20	R1+R2	100
100013260K	CREATININ	JAFFE	2x30+2x30	R1+R2	1.000
101008450K	GAMMA GT	CARBOXI	4x40+1x40	R1+R2	1.250
101294460K	GLUCOSE	GOD-POD	4x60	R1	1.090
101380326K	HDL DIRECT	DIRECT	3x20+1x20	R1+R2	696
101380380K	HDL DIRECT	DIRECT	3x60+1x60	R1+R2	1.040
101016450K	IRON	FERENE	4x40+1x40	R1+R2	1.310
101021450K	LDH	DGKC	4x40+1x40	R1+R2	1.250
101239280K	LDL DIRECT	DIRECT	2X60+1X40	R1+R2	696
101239380K	LDL DIRECT	DIRECT	2x60+1x40	R1+R2	1.040
101319408K	LIPASE	COLOR	40+8	R1+R2	320
100485460K	MAGNESIUM XL	XILIDIL BLUE	4x60	R1	1.090
100384460K	PHOSPHORUS UV	Amm.molibdate	4x60	R1	1.200
100279460K	PROTEIN TOTAL	BIURET	4x60	R1	1.410
100288460K	PROTEIN URYNARY	COOMASSIE	4x60	R1	1.090
101296460K	TRIGLICERIDES	GPO-Trinder	4x60	R1	1.090
101022450K	UREA UV	UREASE	4x40+1x40	R1+R2	870
101012450K	URIC ACID AOX	URICASE	4x40+1x40	R1+R2	1.330
101011460K	URIC ACID MR	URICASE	4x60	R1	1.600



MeDia diagnostici if you need to be sure!

MeDia diagnostici srl, via Costiera 31D/E, 47122 Forlì FC, production unit via Giovanni XXIII, 24, 20090 Milano - Italy
Tel. 0039 0543 782786 fax 0039 0543 782787 email : sales@mediagnostic.biz - www.mediagnostic.biz

Konelab™ 20, T30, T60 packaging with barcode for Konelab™ user defined applications.

Cat.No.	Product description	Method	Kit size	Kit.config.	No. of tests
1011032520K	CHOLINESTERASE	BUTYRRIL	5X20	R1+R2	250
1012924010K	CRP	LEIA	40+10	R1+R2+STD	228
10110393010K	HBA1	TURBI	30+10	R1+R2	160
1019755020K	HOMOCYSTEIN	ENZYMATIC	50+8	R1+R2+std	200
101979220K	HOMOCYSTEIN	ENZYMATIC	25+4	R1+R2+std	200
101969140K	IgA	TURBI	1X40	R1	114
101970140K	IgG	TURBI	1X40	R1	114
101971140K	IgM	TURBI	1X40	R1	114
1012935010H	MICROALBUMINA	LEIA	1X40+1X10	R1+R2	160
1013374010K	RF	LEIA	40+10	R1+R2+STD	160
1013414010K	TRANSFERRIN	TURBI	40+10	R1+R2	160

KEY NON STANDARD ABBREVIATION

AMP	2-amino-2-methyl-1-propanol
BCG	Bromocresol green
BCP	Bromocresol purple
CHOD-PAP	Chol.Oxidase Phenol 4-aminoantipyrine peroxidase
COLOR	Colorimetric
CPC/AMP	O-Cresolphthalein complexone/2-amino-2-methyl-1-propanol
DCA	Dichloroaniline
DEA	Diethanolamine
DPD	Dichlorophenyl diazonium
DGKC	Deutsche Gesellschaft fur Klinische Chemie
ELISA	Enzyme-Linked Immunoassorbent Assay
ECCLS	European Committee for Clinical Laboratory Standards
EPS	Ethylidine PNPG7-blocked Substrate
GLDH	Glutamate Dehydrogenase
GOD-PAP	Glucose Oxidase Phenol 4-Aminoantipyrine Peroxidase
GPO-PAP	Glycerol-3-phosphate oxidase Phenol 4-Aminoantipyrine Peroxidase
IFCC	International Federation for Clinical Chemistry
UV	Ultra violet
UV-GLDH	Ultra Violet- Glutamate Dehydrogenase

REAGENT CHARACTERISTICS

Serum start (one shot)	▶s
Reagent start (one shot)	▶r
All liquid components	L



MeDia diagnostici if you need to be sure !

MeDia diagnostici srl, via Costiera 31D/E , 47122 Forlì FC , production unit via Giovanni XXIII, 24 , 20090 Milano - Italy
Tel. 0039 0543 782786 fax 0039 0543 782787 email : sales@mediagnostic.biz - www.mediagnostic.biz